

*a measurable difference...*



## Cavity Sample Mould

---

### Model: D0005

This four-cavity mould has been produced to fill certain requirements in the specification for Vulcanising Rubber.

### Applications:

- Rubber
- Synthetic Rubber

### Features:

- 4 Cavities
- Hardened stainless steel
- Cavities measuring:  
150mm x 150mm x 2mm deep

### Standards:

- ASTM D3182

### Options:

- Custom made Cavity sizes

### Connections:

- Not Applicable

### Dimensions:

- **H:** 20mm
- **W:** 380mm
- **D:** 380mm
- **Weight:** 10kg



---

10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885 Fax: +61 3 9708 6770

Website: [www.idminstruments.com.au](http://www.idminstruments.com.au) Email: [idm@idminstruments.com.au](mailto:idm@idminstruments.com.au)

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows. Products and details shown may change without notice as IDM Instruments Pty Ltd<sup>®</sup> reserves the right to make alterations to its product range from time to time.

Copyright © 2006 IDM Instruments Pty Ltd<sup>®</sup>. All Rights Reserved.