

*a measurable difference...*



## Circular Cutter

**Model Number: C0032**  
**Cutter, Circular Hand 50cm<sup>2</sup>**

**Model Number: C0038**  
**Cutter, Circular Hand 100cm<sup>2</sup> for Grammage**

The Circular Hand Cutters are fast and easy to operate, creating accurate circular samples in various different mediums. They are available in several different sizes for suitability to many different tests.

### Applications:

- Paper
- Paperboard
- Corrugated Fibreboard
- Solid Fibreboard
- Textiles

### Features:

- Changeable blades
- Easy operation
- 4 pin location points
- Max. Thickness: 7mm
- Safety auto sampling eject plate

### Standards:

- AS • ASTM • AS • BS
- TAPPI • DIN • ISO

### Options:

- Different sized Cutters available upon request

### Connections:

- N/A

### Dimensions:

Varying for different models

- **C0032:** Cutter Flat Crush Test, 50.0cm<sup>2</sup> –0.3cm
- **C0038:** Cutter, Circular Hand 100cm<sup>2</sup> for Grammage



10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885 Fax: +61 3 9708 6770

Website: [www.idminstruments.com.au](http://www.idminstruments.com.au) Email: [idm@idminstruments.com.au](mailto:idm@idminstruments.com.au)

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows. Products and details shown may change without notice as IDM Instruments Pty Ltd<sup>®</sup> reserves the right to make alterations to its product range from time to time.

Copyright © 2006 IDM Instruments Pty Ltd<sup>®</sup>. All Rights Reserved.