

a measurable difference...



Colour Fastness Tester

Model: C0010

This apparatus is one of a series for determining the colour fastness of textiles to various agencies. It determines the colourfastness of textiles to light using an artificial light source.

Applications:

- Fabrics
- Fibres
- Yarn
- NOT polychromatic coloured textiles

Features:

- 500 Watt lamp
- 12 Sample holders
- Black panel thermometer
- Hour Timer
- Countdown Timer

Standards:

- AS 2001.4.21

Options:

- Centering Device

Connections:

- **Electrical:** 220/240 VAC @ 50 HZ or
110 VAC @ 60 HZ
(please specify when ordering)

Dimensions:

- **H:** 400mm • **W:** 650mm • **D:** 500mm
- **Weight:** 10kg



10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885

Fax: +61 3 9708 6770

Website: www.idminstruments.com.au

Email: idm@idminstruments.com.au

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows. Products and details shown may change without notice as IDM Instruments Pty Ltd[®] reserves the right to make alterations to its product range from time to time.

Copyright © 2006 IDM Instruments Pty Ltd[®]. All Rights Reserved.