

Crockmeter

Model: C0019

a measurable difference...

IDM[®]
instruments

The specimen is clamped to the base plate, where a 16mm peg is moved back and forth over the test specimen to determine colour fastness from dry or wet tests.

Applications:

- Textiles
- Carpets

Features:

- Easy manual operation

Benefits:

- Easy to use
- Fast results
- Accurate

Standards:

- BS 1006
- AS 2001.4.3
- AS / NZ 2111.19.1
- AATCC - Test method 8-1972

Options:

- N/A

Connections:

- N/A

Dimensions:

- **H:** 200mm
- **W:** 700mm
- **D:** 220mm
- **Weight:** 20kg

