

*a measurable difference...*



## Cutting Dies - Ruler Type

### Model: C0025

The Ruler Type Cutting Dies are used for cutting of plastic film, paper and rubber samples used for tensile and tear testing (dumbbell, trouser, etc). They are manufactured with a ply backed knife-edge, and can be manufactured to suit various international test standards.

### Application:

- Plastic Film
- Paper
- Rubber

### Features:

- Low volume work

### Standards:

- BS2783
- AS1683.11
- AS2282.7
- ASTM D470
- ASTM D638
- ASTM 1660.2
- ASTM D3574
- AS1301.404
- AS1180.2
- ASTM D412
- ASTM D624
- ASTM D1004
- ASTM D1822

### Options:

- Cutting Press

### Connections:

- N/A

### Dimensions:

- as per requirements



Can be used with the following Cutting Presses



C0022



C0050



S0003

10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885

Fax: +61 3 9708 6770

Website: [www.idminstruments.com.au](http://www.idminstruments.com.au)

Email: [ids@idminstruments.com.au](mailto:ids@idminstruments.com.au)

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows. Products and details shown may change without notice as IDM Instruments Pty Ltd<sup>®</sup> reserves the right to make alterations to its product range from time to time.

Copyright © 2006 IDM Instruments Pty Ltd<sup>®</sup>. All Rights Reserved.