

Cutting Press-2 Tonne Arbour

Model: C0022

a measurable difference...

IDM[®]

instruments

The handle operated manual Arbour Press is used in conjunction with cutting dies to produce testing samples. An optional cutting die adaptor can be used when using the Arbour Press with the Tool Steel cutting dies for fast and repetitious sample preparation.

Applications:

- Paper
- Plastic
- Rubber
- Textiles

Features:

- Ease of operation
- Rack and Pinion
- Table size: 250 mm x 300 mm
- Pressure Foot: 125mm Ø

Benefits:

- Easy to use
- Fast results

Standards:

- N/A

Options:

- Cutting Die Adaptor
- Cutting Dies

Connections:

- N/A

Dimensions:


- **H:** 550mm
- **W:** 350mm
- **D:** 450mm
- **Weight:** 50kg



Can be used with the following Cutting Dies



C0024



C0025



Cutting Die Adaptor