

*a measurable difference...*



## Deep Throat Gauge

---

### Model: T0021

IDM supply a variety of thickness measuring devices for various applications, one of which being the Deep Throat Gauge.

The Deep Throat Gauge has been specifically designed for measuring the thickness at a greater distance from the samples edge.

### Application:

- Plastics
- Textiles
- Rubber

### Specifications:

- Depth: 450mm
- Open to: 10mm
- Read to: 0.01mm

### Standards:

- N/A

### Options:

- Digital or Analogue gauge
- Anvil sizes
- Throat sizes
- Statistical Printer (with Digital Gauge option)

### Connections:

- N/A

### Dimensions:

- **H:** 210mm
- **W:** 600mm
- **D:** 32mm
- **Weight:** 5kg



---

10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885

Fax: +61 3 9708 6770

Website: [www.idminstruments.com.au](http://www.idminstruments.com.au)

Email: [idm@idminstruments.com.au](mailto:idm@idminstruments.com.au)

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows. Products and details shown may change without notice as IDM Instruments Pty Ltd<sup>®</sup> reserves the right to make alterations to its product range from time to time.

Copyright © 2006 IDM Instruments Pty Ltd<sup>®</sup>. All Rights Reserved.