

*a measurable difference...*



## Dry Bath Aging Block

---

### Model: D0001

To be used for determining the long-term effects of elevated temperature on various materials like rubber, plastics, textile etc. Samples are placed in glass tubes. The 24 holes for insertion of the tubes are 38mm by 280mm.

### Applications:

- Rubber
- Plastic
- Textiles

### Features:

- Set of glassware
- Temperature controller: Ambient – 300°C

### Standards:

- Not Applicable

### Options:

- Not Applicable

### Connections:

- **Electrical:** 220/240 VAC @ 50 HZ or 110 VAC @ 60 HZ

*(please specify when ordering)*

### Dimensions:

- **H:** 500mm
- **W (OD):** 565mm
- **Control Box:** 150mm
- **Weight:** 90kg



---

10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885

Fax: +61 3 9708 6770

Website: [www.idminstruments.com.au](http://www.idminstruments.com.au)

Email: [idm@idminstruments.com.au](mailto:idm@idminstruments.com.au)

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows. Products and details shown may change without notice as IDM Instruments Pty Ltd<sup>®</sup> reserves the right to make alterations to its product range from time to time.

Copyright © 2006 IDM Instruments Pty Ltd<sup>®</sup>. All Rights Reserved.