

a measurable difference...



Drying Cabinet

Model: D0012

The Drying Cabinet is used for the rapid drying of textile samples. The cabinet is fitted with three rows of folding pull-out hangers, corresponding to a total of 16 metres of clothesline to support the samples on. The Drying Cabinet can be operated at different temperatures, ranging from Normal, Medium or Air Fluff temperature.

Application:

- Drying of Textile Samples

Features:

- Racks equivalent to 16m
- Racks also fitted inside of door
- Dehumidification capacity: 14 g/min
- Airing capacity: 180 m³/h

Benefits:

- The Drying Cabinet has multiple temperature ranges allowing you to simulate both “normal” drying regimes or rapidly dry items for maximum test efficiency
- Our clever system racking creates 16 meters of effective drying space, so you can dry more items in this small footprint cabinet
- 8 preset timers offering “set and forget” ease of use

Standards:

- N/A

Options:

- Stainless Steel Construction

Connections:

- **Energy consumption:** Normal 0.9 kWh/kg

- **Electrical:** 220/240 VAC @ 50 HZ or
110 VAC @ 60 HZ

(please specify when ordering)

Dimensions:

- **H:** 1,730 mm
- **W:** 595 mm
- **D:** 640 mm
- **Weight:** 60 kg



10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885 Fax: +61 3 9708 6770

Website: www.idminstruments.com.au Email: idm@idminstruments.com.au

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows. Products and details shown may change without notice as IDM Instruments Pty Ltd[®] reserves the right to make alterations to its product range from time to time.