

Manual Sample Cutters

Model: C0028, C0029, C0031

a measurable difference...

IDM[®]
instruments

Ring Crusher Sample Cutter: C0028

STFI/ Tensile: C0029

Tear Test Sample Cutter: C0031

Precision Sample Cutters for cutting testing samples for paper to the following sizes as standards.

IDM can also manufacture custom sizes if required.

Applications:

- Paper

Features:

- C0028: 12.7 x 152 mm
- C0029: 15 x 300 mm
- C0031: 50 x 62 mm
- Can be resharpened

Benefits:

- Easy to use
- Fast results
- Accurate

Standards:

- C0028: AS 1301.407
- C0029: AS 1301.448
- C0031: AS 1301.400

Options:

- TBA

Connections:

- N/A

Dimensions:

- As per requirements

