

*a measurable difference...*



## Manual Sample Cutters

---

### Models:

**Ring Crush Sample Cutter:** C0028  
**STFI / Tensile:** C0029  
**Tear Test Sample Cutter:** C0031

Precision Sample Cutters for cutting testing samples for paper to the following sizes as standards. IDM can also manufacture custom sizes if required.

### Application:

- Paper

### Features:

- C0028: 12.7 x 152 mm
- C0029: 15 x 300 mm
- C0031: 50 x 62 mm
- Can be resharpened

### Standards:

- C0028: AS 1301.407
- C0029: AS 1301.448
- C0031: AS 1301.400

### Options:


- TBA

### Connections:

- N/A

### Dimensions:

- As per Cutter requirements



---

10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885 Fax: +61 3 9708 6770

Website: [www.idminstruments.com.au](http://www.idminstruments.com.au) Email: [idm@idminstruments.com.au](mailto:idm@idminstruments.com.au)

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows. Products and details shown may change without notice as IDM Instruments Pty Ltd<sup>®</sup> reserves the right to make alterations to its product range from time to time.

Copyright © 2006 IDM Instruments Pty Ltd<sup>®</sup>. All Rights Reserved.