

a measurable difference...



Pneumatic Sample Cutter

Model: C0043

The Pneumatic Sample Cutter is a fast and easy to operate mechanism used for sample preparation.

The device is intended for preparation of circular / square / rectangle test pieces with uniform edges. It is pneumatic, for easy operation, and repetitious sample cutting.

Applications:

- Paper
- Cartonboard

Features:

- Front paper load for easy operation
- One or more sheets can be cut simultaneously
- Punch and Die made from tempered alloy steel
- Safety Guard
- Grammage range: 35 to 1,000gsm
- Thickness: Up to 4mm thick

Standards:

- AS
- TAPPI
- ISO
- BS
- SCAN
- APPITA
- DIN

Options:

- Cutter sizes selection:
 - 100cm²
 - 38.1x70mm
 - 62x50mm
 - 63x53mm
 - 50x62mm
 - 38x70mm
 - 15x100mm
- Custom cutter size available upon request.
- Foot Switch



Connections:

- Air supply: 80 — 100 psi

Dimensions:

- H: 250mm
- W: 230mm
- D: 380mm
- Weight: 23kg

10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885

Fax: +61 3 9708 6770

Website: www.idminstruments.com.au

Email: idm@idminstruments.com.au

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows. Products and details shown may change without notice as IDM Instruments Pty Ltd[®] reserves the right to make alterations to its product range from time to time.

Copyright © 2006 IDM Instruments Pty Ltd[®]. All Rights Reserved.