

a measurable difference...



Polarising Box

Model: P0008

The Polarising Box is used to view the stress points created during manufacture of clear plastic and glass products.

The samples can also be extruded film. The Polarising Box has two open sides so long items can be passed through.

Application:

- Plastic films
- Glass

Features:

- Polarising Filter mounted at the top for viewing
- Light Source
- Fluorescent Tubes, mounted at base with a diffuser as the sample support panel

Standards:

- N/A

Options:

- Larger Polarising Box (Consult Factory)

Connections:

- **Electrical:** 220/240 VAC @ 50 HZ or 110 VAC @ 60 HZ

(please specify when ordering)

Dimensions:

- **H:** 400mm
- **W:** 400mm
- **D:** 400mm
- **Weight:** 7kg



10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885

Fax: +61 3 9708 6770

Website: www.idminstruments.com.au

Email: idm@idminstruments.com.au

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows. Products and details shown may change without notice as IDM Instruments Pty Ltd[®] reserves the right to make alterations to its product range from time to time.

Copyright © 2006 IDM Instruments Pty Ltd[®]. All Rights Reserved.