

# Sample Cutting Press-10 Tonne Hydraulic

*a measurable difference...*

# IDM<sup>®</sup>

instruments

Model: S0003

The Sample Cutting Press 10 Tonne Hydraulic has been designed to be used for sample preparation for cutting of rubber, plastic and carton board. It has been manufactured with two spring-return safety switches, which need to be pushed and held to activate the cutting action. Tool Steel or Ruler Type Cutting Dies can be used with the Cutting Press for fast and repetitious cutting.

## Applications:

- Sample cutting of Rubber, Plastic and Carton Board

## Features:

- Electric - Hydraulic pump
- Cutting Stroke: 125mm
- Throat: 125mm
- Cutting Base: 300 x 300mm

## Benefits:

- Easy to use
- Fast results
- Accurate

## Standards:

- N/A

## Options:

- Cutting Dies to be ordered separately  
(please refer to C0024 & C0025 data sheets)
- Cutting Die Adaptor

## Connections:

- **Electrical:** 220/240 VAC @ 50 HZ or  
110 VAC @ 60 HZ  
(please specify when ordering)

## Dimensions:

- **H:** 300mm
- **W:** 350mm
- **D:** 600mm
- **Weight:** 60kg



Cutting Die Adaptor



Can be used with the following Cutting Dies



C0024



C0025