

*a measurable difference...*



## Sample Cutting Press 10 Tonne Hydraulic

**Model: S0003**

The Sample Cutting Press 10 Tonne Hydraulic has been designed to be used for sample preparation for cutting of rubber, plastic and carton board. It has been manufactured with two spring-return safety switches, which need to be pushed and held to activate the cutting action. Tool Steel or Ruler Type Cutting Dies can be used with the Cutting Press for fast and repetitious cutting.

### Application:

- Sample cutting of Rubber, Plastic and Carton Board

### Features:

- Electric - Hydraulic pump
- Cutting Stroke: 125mm
- Throat: 125mm
- Cutting Base: 300 x 300mm

### Standards:

- N/A

### Options:

- Cutting Dies to be ordered separately (please refer to C0024 & C0025 data sheets)
- Cutting Die Adaptor

### Connections:

- **Electrical:** 220/240 VAC @ 50 HZ or 110 VAC @ 60 HZ (please specify when ordering)

### Dimensions:

- **H:** 300mm    • **W:** 350mm    • **D:** 600mm
- **Weight:** 60kg




Cutting Die Adaptor

Can be used with the following Cutting Dies



C0024



C0025

10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885

Fax: +61 3 9708 6770

Website: [www.idminstruments.com.au](http://www.idminstruments.com.au)

Email: [idm@idminstruments.com.au](mailto:idm@idminstruments.com.au)

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows. Products and details shown may change without notice as IDM Instruments Pty Ltd<sup>®</sup> reserves the right to make alterations to its product range from time to time.