

a measurable difference...



Schopper Riegler Pneumatic

Model: SR/P

The SR/P is used for the determination of drainability of pulp suspension in water.

The drainability is the measurable index of the refining degree of pulps. The refining of pulps is one of the most important stages in the paper production process and influences strongly the sheet forming and its physical properties.

Applicable to all kinds of pulp in aqueous suspension.



SR-DR

MAJOR TECHNICAL CHARACTERISTICS	
Speed of the sealing cone	(100 ± 10) mm/s
Calibrated bottom orifice	(149,0 ± 1,0) s for 1 L of H ₂ O
Drainage chamber volume	1000 mL
Drainage area	100 cm ²
Pneumatic drive	2,0 kgf/cm ²
Measuring scales	(0 to 100)*SR
Resolution	1°SR

The SR/P also features iron-cast base; large precision circular level; leveling feet; stainless steel structure; chrome plated bronze drainage chamber, funnel, spreader cone and sealing cone; plexiglas safeguard; handle for drainage chamber; ISO standard pneumatic components; hook spanner and fixing base for comfortable exchange of phosphor bronze wire screen.

Optional accessories:

· Filter for water model GS200UV - 5 L/min according to ISO 14487

Applicable standards:

· ISO 5267/1

Air supply	2 kgf/cm ² instrument quality
Air consumption	0,1 m ³ /h
Dimensions	(330 x 470 x 1440) mm (W x L x H)
Weight	25 kg

Note: due to constant development our equipment design and specifications are subject to change without notice.

10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885

Fax: +61 3 9708 6770

Website: www.idminstruments.com.au

Email: ids@idminstruments.com.au

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows. Products and details shown may change without notice as IDM Instruments Pty Ltd[®] reserves the right to make alterations to its product range from time to time.