

An edition of our latest news...

EXHIBITION NOTICES:**Appita 2007: 61st Appita Annual Conference & Exhibition**

Venue: Gold Coast Convention & Exhibition Centre, Queensland, AUSTRALIA
Dates: 6th May – 9th May 2007
Booth # 35

Furnitex 2007

Venue: Melbourne Exhibition Centre, Victoria, AUSTRALIA
Dates: 19th July – 22nd July 2007
Booth # M46

NEW DEVELOPMENT:**Fabric Vertical Burn Tester**

The Fabric Vertical Burn Tester was manufactured to test the burning characteristics of vertically oriented fabrics, a useful quality control tool used in the laboratory. The added benefit of the IDM Fabric Vertical Burn Tester is the ability to modify and manufacture the instrument to be able to conform to various international standards. Rather than using 2 or 3 burner trays that allows different burner angle positions for ignition, IDM have manufactured 1 burner tray with an adjustable burner, allowing for easy angle positioning of the burner to accommodate different burning testing positions. For added health and safety of operators, IDM have also manufactured a Stainless Steel Enclosure for the Fabric Vertical Burn Tester, enabling the test to be operated inside the enclosure at all times.

Thickness Gauge Digital Indicator & Stand

The Thickness Gauge Digital Indicator & Stand can be used on many different applications where an accurate measurement of thickness is required. It can be modified to be used with different type digital gauges offering outputs in mm/inch measurements, and different ranges available. A dial gauge can also be fitted to the Thickness Gauge where applicable. Added benefits of the Thickness Gauge is the ability to output results onto an external printer while using a thumb trigger for easy testing applications.

February 2007

Website:

www.idminstruments.com.au

Email:

idm@idminstruments.com.au

10-11 Colrado Court, Hallam, VIC 3803, Australia
Phone: 61 3 9708 6885 Fax: 61 3 9708 6770

Page 1