

# Deep Throat Gauge

Model: T0021

*a measurable difference...*

**IDM**<sup>®</sup>  
instruments

IDM supply a variety of thickness measuring devices for various applications, one of which being the Deep Throat Gauge.

The Deep Throat Gauge has been specifically designed for measuring the thickness at a greater distance from the samples edge.

## Applications:

- Plastics
- Textiles
- Rubber

## Specifications:

- Depth: 450mm
- Open to: 10mm
- Read to: 0.01mm

## Benefits:

- Easy to use
- Fast results
- Accurate

## Standards:

- N/A

## Options:

- Digital or Analogue gauge
- Anvil sizes
- Throat sizes
- Statistical Printer (with Digital Gauge option)

## Connections:

- N/A

## Dimensions:

- **H:** 210mm
- **W:** 600mm
- **D:** 32mm
- **Weight:** 5kg

